Programat® CS4

Valid from software V2.1

Program structure

Program	Name	Description	Duration[6]	Programat® Sinter Tray	IPS e.max® CAD Crystallization Tray	IPS Speed Tray
1	IPS e.max ZirCAD Prime Speed sintering	Program for the fast sintering of IPS e.max ZirCAD Prime[1], incl. pre-drying	27 min[2]	✓[3]	–	–
2	IPS e.max ZirCAD MT Multi Speed sintering crown	Program for the fast sintering of IPS e.max ZirCAD MT Multi[1] crowns, incl. pre-drying	60 min[2]	✓[3]	–	–
3	IPS e.max ZirCAD MT Multi Sintering crown/bridge	Program for sintering of IPS e.max ZirCAD MT Multi[1] crowns/bridges with higher translucency, incl. pre-drying	114 min[2]	✓[3]	–	–
4	IPS e.max ZirCAD LT Speed sintering	Program for the fast sintering of IPS e.max ZirCAD LT[1], incl. pre-drying	27 min[2]	✓[3]	–	–
5	IPS e.max ZirCAD LT Sintering	Program for sintering of IPS e.max ZirCAD LT[1] with higher translucency, incl. pre-drying	47 min[2]	✓[3]	–	–
6	IPS e.max ZirCAD Stains/Glaze firing	Program for the Stains/ Glaze firing of IPS e.max ZirCAD with IPS e.max CAD Crystall./ materials	19 min	✓[3]	✓	–
7	IPS e.max CAD Speed crystallization (HT, MT, LT)	Program for the fast crystallization of IPS e.max CAD HT, MT or LT	17 min	✓[3]	–	✓
8	IPS e.max CAD Crystallization (HT, MT, LT)	Program for the crystallization of IPS e.max CAD HT, MT or LT	26 min	✓[3]	✓	–
9	IPS e.max CAD Crystallization (MO, Impulse)	Program for the crystallization of IPS e.max CAD MO or Impulse	33 min	✓[3]	✓	–
10	IPS e.max CAD Corrective/Stains/Glaze firing	Program for the Corrective, Stains and Glaze firing of IPS e.max CAD with IPS e.max CAD Crystall./ materials	22 min	✓[3]	✓[4]	✓[4]
11	IPS Empress CAD Stains/Glaze firing	Program for the Stains and Glaze firing of IPS Empress CAD with IPS e.max CAD Crystall./ materials	22 min	✓[3]	✓	–
12 –30	Individual programs			✓[3]	(✓)[5]	(✓)[5]

[1] This sinter program may only be used for restorations made from blocks.
[2] Duration without pre-drying. Duration is prolonged by 8–10 min when including pre-drying time.
[3] Programat Sinter Tray must always be used
[4] Use IPS e.max CAD Crystallization Tray or IPS Speed Tray
[5] Use a crystallization tray suitable for the program type. Not suitable for sintering programs.
[6] R&D Ivoclar, Schaan

Important: Predrying is mandatory for crystallization and glaze programs.

For more detailed information, please refer to the Instructions for Use of the respective material.

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