Programat® CS6
Software Revision Log

Version	Changes from previous version	Status
V2.2 26.01.2026	**Software optimized** • Vacuum test optimized • Threshold values for error messages optimized • Valves switch during self-test • Various software optimizations	Currently delivered version
V2.1 12.09.2023	**Software optimized** • Programs for "IPS e.max ZirCAD Prime Block" added • Program order modified **Ivoclar programs will be overwritten!** • Not possible pre-drying options deactivated • Beep at the end of programs • Text of the cooling screen adjusted • Various software optimizations	No longer available
V2.0 26.07.2022	**Software optimized** • Waiting time between two programs was shortened • New program can be started immediately after the previous program is finished • New Ivoclar-Logo implemented • Various software optimizations	No longer available
V1.3 21.01.2022	**Software optimized** • Parameters of program P1 and pre-drying "cryst." optimized	No longer available
V1.2 06.09.2021	**Software optimized** • Software-bug was fixed. Ivoclar programs will be set back to factory values. Individual programs are not affected	No longer available
V1.1 30.06.2021	**Basic version**	No longer available

Version 26.01.2026